

IHE Integration Statement			
Vendor	Product Name	Version	Date
CGM Clinical Österreich	Integration Package (i.Pack)	4.1.x	6/2/17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Application	None
Consistent Time	Time Client	None
Pharmacy Dispense	Content Consumer	View Option
Pharmacy Dispense	Content Consumer	None
Document Metadata Subscription	Document Metadata Notification Recipient	None
Document Metadata Subscription	Document Metadata Subscriber	None
Document Metadata Subscription	Document Metadata Subscriber	Patient-Independent Subscription
eCR	eCR Client	None
Pharmacy Pharmaceutical Advice	Content Consumer	View Option
Pharmacy Pharmaceutical Advice	Content Consumer	None
Patient Administration Management	Patient Demographics Supplier	None
Patient Administration Management	Patient Demographics Supplier	Link / Unlink
Patient Administration Management	Patient Demographics Supplier	Merge
Patient Administration Management	Patient Encounter Supplier	None
Patient Demographics Query	Patient Demographics Consumer	None

Internet address for vendor's IHE information :

http://www.cgm.com/at/krankenhaus_klinik_1/kommunikationsloesungen_1/cgm_ihe_1/downloads_8/downloads_ihe.de.jsp

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
CGM Clinical Österreich	Integration Package (i.Pack)	4.1.x	6/2/17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identity Source	None
Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-Reference HL7 V3	Patient Identity Source	None
Pharmacy Prescription	Content Consumer	View Option
Pharmacy Prescription	Content Consumer	None
Cross-Enterprise Document Reliable Interchange	Document Source	None
Cross-Enterprise Document Reliable Interchange	Metadata-Limited Document Source	None
Cross-Enterprise Document Sharing	Document Administrator	None
Cross-Enterprise Document Sharing	Document Administrator	Update Metadata

Internet address for vendor's IHE information :

http://www.cgm.com/at/krankenhaus_klinik_1/kommunikationsloesungen_1/cgm_ihe_1/downloads_8/downloads_ihe.de.jsp

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
CGM Clinical Österreich	Integration Package (i.Pack)	4.1.x	6/2/17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing	Document Administrator	DEPRECATED>Delete Metadata
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Source	None
Cross-Enterprise Document Sharing	Document Source	Document Replacement
Cross-Enterprise Document Sharing	Document Source	Document Transformation
Cross-Enterprise Document Sharing	Document Source	Document Addendum
Cross-Enterprise Document Sharing	Document Source	Folder Management
Cross-Enterprise Document Sharing for Imaging	Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	View Option
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	None
Cross Enterprise User Assertion	X-Service User	None
Cross Enterprise User Assertion	X-Service User	Subject-Role

Internet address for vendor's IHE information :

http://www.cgm.com/at/krankenhaus_klinik_1/kommunikationsloesungen_1/cgm_ihe_1/downloads_8/downloads_ihe.de.jsp

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
CGM Clinical Österreich	Integration Package (i.Pack)	4.1.x	6/2/17

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross Enterprise User Assertion	X-Service User	PurposeOfUse

Internet address for vendor's IHE information :

http://www.cgm.com/at/krankenhaus_klinik_1/kommunikationsloesungen_1/cgm_ihe_1/downloads_8/downloads_ihe.de.jsp

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net