

Control #:	Control Name					
Test	Acc #	Run Date	Result	Flag	Corrective Action	Run By
Chemistry						
.C1 Low Chemistry						
Glucose		Units: mg/dL				
	1175	1/18/2013 9:27 AM	67		Result out of range. Retested sample and result was within limits. No action required, SL	LAB
	3753	3/20/2013 9:32 AM	68		WEEKLY REVIEW OF QC	LAB
BUN		Units: mg/dL				
	1186	1/20/2013 9:47 AM	5		recalibrated inst	LAB
	3744	3/19/2013 9:35 AM	10	Outside 3SD	Reran	LAB
Sodium		Units: mEq/L				
	1149	1/12/2013 9:39 AM	124	Two Consecutive Outside 2SD	Result out of range; retested and result was within expected limits. No further action at this time.	LAB
	3661	3/3/2013 10:15 AM	128		recalibrated-br	LAB
	3664	3/4/2013 10:08 AM	130	Four Consecutive Outside 1SD	reran QC- br	LAB
	3753	3/20/2013 9:32 AM	128		WEEKLY REVIEW	LAB
Potassium		Units: mEq/L				
	1149	1/12/2013 9:39 AM	2.6	Four Consecutive Outside 1SD	Result was out of range. Repeat test was within expected limits. No further action at this time.	LAB
	1197	1/24/2013 10:04 AM	2.6	Four Consecutive Outside 1SD	Result out of range. Retested sample and result was within expected limits. No further action required.	LAB
Hematology						
.H2 Hematology Abnormal						
WBC		Units: K/ μ L				
	1196	3/5/2013 9:53 AM	3.0	Outside 3SD	Reran controls CS	LAB