



IHE Integration Statement			
Vendor	Product Name	Version	Date
Systema Human Information Systems	Integration Package (i.Pack)		7/14/14

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Node	None
Consistent Time	Time Client	None
Patient Demographics Query	Patient Demographics Consumer	None
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identity Source	None
Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-Reference HL7 V3	Patient Identity Source	None
Cross-Enterprise Document Sharing	Document Administrator	None
Cross-Enterprise Document Sharing	Document Administrator	Update Metadata
Cross-Enterprise Document Sharing	Document Administrator	Delete Metadata
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Source	None
Cross-Enterprise Document Sharing	Document Source	Folder Management
Cross-Enterprise Document Sharing	Document Source	Document Addendum
Cross-Enterprise Document Sharing	Document Source	Document Transformation
Cross-Enterprise Document Sharing	Document Source	Document Replace

Internet address for vendors IHE information : http://www.systema.info	
Links to Standards Conformance Statements for the Implementation	
HL7	
DICOM	
Link to general information on IHE : www.ihe.net	



IHE Integration Statement			
Vendor	Product Name	Version	Date
Systema Human Information Systems	Integration Package (i.Pack)		7/14/14

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing for Imaging	Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None
Cross Enterprise User Assertion	X-Service User	None
Cross Enterprise User Assertion	X-Service User	PurposeOfUse
Cross Enterprise User Assertion	X-Service User	Subject-Role

Internet address for vendors IHE information : http://www.systema.info	
Links to Standards Conformance Statements for the Implementation	
HL7	
DICOM	
Link to general information on IHE : www.ihe.net	