

FIREWALLANFORDERUNGEN

VPN-ZUGANGSDIENST CGM

Schnittstelle	URL	IPs / IP-Ranges	Port/s
DNS Server	ns1.telemed-ti.net ns2.telemed-ti.net ns3.telemed-ti.net	185.39.252.249 185.39.253.249 176.28.52.198	53/udp 53/tcp
VPN Konzentratoren	*.f-vpnzugd.telemed-ti.net *.d-vpnzugd.telemed-ti.net	185.188.0.11 - 185.188.0.46 185.188.2.11 - 185.188.2.46	esp 500/udp 4500/udp
SIS	vpnk-si-01.f-vpnzugd.telemed-ti.net vpnk-si-02.f-vpnzugd.telemed-ti.net vpnk-si-01.d-vpnzugd.telemed-ti.net vpnk-si-02.d-vpnzugd.telemed-ti.net	185.188.0.71 185.188.0.72 185.188.2.71 185.188.2.72	esp 500/udp 4500/udp
Registrierungsserver	https://register.f-vpnzugd.telemed-ti.net:8443/RegistrationServer/services/provisioningPort https://register.d-vpnzugd.telemed-ti.net:8443/RegistrationServer/services/provisioningPort	185.188.0.4 185.188.2.4	8443/tcp
Hash&URL Server	https://hashandurl.f-vpnzugd.telemed-ti.net/ https://hashandurl.d-vpnzugd.telemed-ti.net/	185.188.0.7 185.188.2.7	80/tcp