

IHE Integration Statement			
Vendor	Product Name	Version	Date
CGM LAB Belgium	MOLIS	4.41	12/08/20

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Application	ATX: FHIR FEED
Audit Trail and Node Authentication	Secure Application	None
Audit Trail and Node Authentication	Secure Application	STX: No Secure Transport
Consistent Time	Time Client	None
Laboratory Code Set Distribution	Code Set Master	Batch Options
Laboratory Code Set Distribution	Code Set Master	None
Patient Administration Management	Patient Demographics Consumer	Merge
Patient Administration Management	Patient Demographics Consumer	None
Patient Administration Management	Patient Encounter Consumer	None
Sharing Lab Report	Content Creator	None
Cross-Enterprise Document Media Interchange	Portable Media Creator	USB Document Storage
Cross-Enterprise Document Media Interchange	Portable Media Creator	None
Cross-Enterprise Document Media Interchange	Portable Media Creator	ZIP file sent over email