

**CGM TELEMEDICINE** is the division of CompuGroup Medical group specialized in the development of solutions and technologies in the telemedicine field with the aim of simplifying the relationship with healthcare professionals, improving accessibility and quality of remote assistance.

**PRODUCER**

CGM TELEMEDICINE



**CONTACTS**

 [commerciale.buhs.it@cgm.com](mailto:commerciale.buhs.it@cgm.com)

**WEB SITE**




Synchronizing Healthcare



**CompuGroup  
Medical**

# CGM HI 3 LEADS ECG



**CGM | TELEMEDICINE**  
Connected Patient Solutions



**CompuGroup  
Medical**


# TELECARDIOLOGY

From diagnosis to home telemonitoring, with **CGM HI 3 LEADS ECG** it is possible wherever you are.

A single device, with three leads, for different ways of use.



- Lightweight and comfortable to wear
- Patient app
- Digital diary
- Activity's levels
- Class IIa medical devices
- Standard patches
- IP67
- Battery life: up to 7 days



By using Bluetooth Connection, **CGM HI 3 LEADS ECG** sent ECG signals to the patient app.

ECG traces will be sent to the web platform and will be available for visualization and analysis to the clinicians

## How to use



### TELE HOLTER

Remote Holter up to 30 days with web-based platform data transmission and dedicated advanced holter analysis software



### TELEMONITORING

Remote monitoring up to 7 days with online transmission and data analysis in a web-based platform by using an integrated and dedicated algorithm



### TELE ECG

Remote rest ECG examination when symptoms occur or on scheduled sessions, with online transmission and data analysis in web-based platform

**CGM HI 3 LEADS ECG** is Integrated with the platform **CGM CARE MAP**, web solution that allows healthcare personnel to always keep in touch with their patients.

Thanks to this integration, you will have the possibility to monitor not only cardiac rhythm but also several vital parameters (oxygen saturation, blood pressure) which help to have a complete clinical picture of your patient's health.

**CGM HI 3 LEADS ECG** revolutionizes and simplifies dispensing methods of telecardiology services and adapts to several application fields:

- Diagnosis main cardiac arrhythmias
- Cardiological evaluation during physical activity and rehabilitation
- Therapeutical response's evaluation
- Pre and postsurgery telemonitoring
- Chronic patient monitoring
- Post discharge telemonitoring