

IHE Integration Statement			
Vendor	Product Name	Version	Date
CompuGroup Medical Software GmbH (formerly CGM Clinical Österreich)	Integration Package (i.Pack)_2025	5.0.2	2/3/26

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Application	None
Audit Trail and Node Authentication	Secure Application	STX:TLS 1.2 Floor using BCP195
Audit Trail and Node Authentication	Secure Application	ATX: TLS Syslog
Consistent Time	Time Client	None
Document Metadata Subscription	Document Metadata Notification Recipient	None
Document Metadata Subscription	Document Metadata Subscriber	None
Document Metadata Subscription	Document Metadata Subscriber	Patient-Independent Subscription
Mobile access to Health Documents	Document Recipient	None
Mobile access to Health Documents	Document Recipient	XDS on FHIR
Mobile access to Health Documents	Document Responder	None
Mobile access to Health Documents	Document Responder	XDS on FHIR
Patient Demographics Query	Patient Demographics Consumer	None
Patient Demographics Query for Mobile	Patient Demographics Consumer	None
Patient Demographics Query for Mobile	Patient Demographics Consumer	Patient Search

Internet address for vendor's IHE information :
https://www.cgm.com/aut_de/produkte/intersektoral/cgm-ihe.html

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
CompuGroup Medical Software GmbH (formerly CGM Clinical Österreich)	Integration Package (i.Pack)_2025	5.0.2	2/3/26

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None
Patient Identifier Cross-referencing	Patient Identity Source	None
Patient Identifier Cross-referencing	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-referencing for mobile	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-referencing HL7 V3	Patient Identity Source	None
Patient Identifier Cross-referencing HL7 V3	Patient Identifier Cross-reference Consumer	None
Web-based Image Access	Imaging Document Consumer	None
Cross-Enterprise Document Reliable Interchange	Document Source	None
Cross-Enterprise Document Reliable Interchange	Metadata-Limited Document Source	None
Cross-Enterprise Document Sharing	Document Administrator	None
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Source	Folder Management

Internet address for vendor's IHE information :
https://www.cgm.com/aut_de/produkte/intersektoral/cgm-ihe.html

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
CompuGroup Medical Software GmbH (formerly CGM Clinical Österreich)	Integration Package (i.Pack)_2025	5.0.2	2/3/26

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing	Document Source	None
Cross-Enterprise Document Sharing	Document Source	Document Addendum
Cross-Enterprise Document Sharing	Document Source	Document Transformation
Cross-Enterprise Document Sharing	Document Source	Document Replacement
Cross-Enterprise Document Sharing	Integrated Document Source/Repository	None
Cross-Enterprise Document Sharing	Integrated Document Source/Repository	Document Replacement
Cross-Enterprise Document Sharing for Imaging	Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	None
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	View Option
Cross Enterprise User Assertion	X-Service User	Subject-Role
Cross Enterprise User Assertion	X-Service User	None

Internet address for vendor's IHE information :
https://www.cgm.com/aut_de/produkte/intersektoral/cgm-ihe.html

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
CompuGroup Medical Software GmbH (formerly CGM Clinical Österreich)	Integration Package (i.Pack)_2025	5.0.2	2/3/26

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross Enterprise User Assertion	X-Service User	PurposeOfUse

Internet address for vendor's IHE information :
https://www.cgm.com/aut_de/produkte/intersektoral/cgm-ihe.html

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net