

## Dedicated Scanning Server—6 HDD (2U)



**Rackmount** 2U

**Processor** Intel® XEON® Silver 4216—Hexadeca Core

**Memory** 32GB RAM

**Storage** RAID-1 mirrored pairs (2 - 4TB, 2 - 2TB, 2 - 1TB)

**Operating System** Windows® Server 2019 Standard—16 Core

Time to upgrade?  
Microsoft ends extended support of Windows® 8 and Server 2012 in 2023.

**Ideal Use** Large labs or any lab using the CGM LABDAQ® Scanning/Attachments option

## Dedicated Server—6 HDD (2U)



**Rackmount** 2U

**Processor** Intel® XEON® Silver 4216—Hexadeca Core

**Memory** 32GB RAM

**Storage** RAID-1 mirrored pairs (6 - 1TB)

**Operating System** Windows® Server 2019 Standard—16 Core

**Ideal Use** Large labs with 3 or more remote lab locations, or single labs with >15 workstations and >3 file-based analyzers

## Dedicated Server—2 HDD\*



**Processor** Intel® XEON® E-2136 v6—Hexa Core

**Memory** 32GB DDR4 RAM

**Storage** Mirrored (2 - 1TB SSD) & 250GB NVMe M.2 Boot Drive

**Operating System** Windows® Server 2019 Standard—16 Core

*\*Optional Rackmount 1U with Mirrored (2 - 1TB HDD; hot swappable)*

**Ideal Use** Mid-sized labs with 5-15 workstations; any lab with 2-3 file-based analyzers; facilities with 1-2 remote labs and only 1 unidirectional analyzer connection per remote lab

## Server-Workstation



**Processor** AMD Ryzen™ 5 3400G—Quad Core

**Memory** 32GB DDR4 RAM

**Storage** 1TB SSD & 250GB NVMe M.2 Boot Drive  
External USB Drive

**Operating System** Windows® 10 Pro

## Workstation

**Processor** AMD Ryzen™ 5 3400G—Quad Core

**Memory** 16GB DDR4 RAM

**Storage** 250GB NVMe M.2 Drive

**Operating System** Windows® 10 Pro

**Ideal Use** Small labs with 1-4 workstations, 1 file-based analyzer, 1 EMR interface (orders and results), and 1-2 reference lab interfaces



The lab division at CompuGroup Medical makes setting up your laboratory easy. Choose pre-configured hardware designed to meet or exceed the system requirements of **CGM LABDAQ®** and the specific needs of your lab. Choose the appropriate hardware, and enjoy the CGM advantage. For more information regarding system requirements, visit [www.cgmus.com/sys\\_reqs/](http://www.cgmus.com/sys_reqs/).

With CGM hardware, your lab will enjoy:

- Hardware designed to system requirements
- First year hardware warranty at no charge
- First year hardware support at no charge
- Server migration
- Technical support assistance when reconnecting analyzers to workstations
- In-stock availability

## LABEL PRINTERS



ZD410 Ethernet  
Network

**Ideal Use**  
Multi-site  
facilities, or  
labs with  
three or more  
workstations



ZD410 USB  
Non-Network

**Ideal Use**  
Single-site  
facilities with  
one or two  
workstations



**Contact us  
to learn more**

1.800.359.0911 opt. 3

[sales.lab.us@cgm.com](mailto:sales.lab.us@cgm.com)  
[cgm.com/labdaq](http://cgm.com/labdaq)