



Dear Valued Customer,

I am pleased to announce the General Availability of CGM LABDAQ® version 20.7 from CompuGroup Medical.

The key improvements with this release that have been introduced since version 20.5 include:

For all laboratories:

Inventory Improvements

- A new inventory type of 'Equipment' has been added to Inventory Type, for more robust inventory tracking
- An attachments feature has been added to Inventory Maintenance to help track any associated attachments of inventory items
- An inventory group filter has been added to help more efficiently review and track inventory
- Consumable and Equipment can now be added from the '*Workflow Management*' screen (see also below)

Reporting Improvements

- In Batch Report Printing, 'Accession', 'Received Date' and 'Collection Date' have been added as sort options
- In Batch Label Printing, 'Sort Order' has been added
- Corrected result history is now viewable on the Result Report for all prior corrected results for a specific test
- With memorized report sharing, users in a group can share memorized reports with other User Groups and/or Users

Controls Improvements

- Control Tests now follow Test Sorting Setup (rather than being alphabetized)

Requisition Improvements

- Improved container management, where containers can now be received through the 'Sample Receiving' window, and tracked via the 'Received' tab of Requisition Maintenance
- Improved usability and tracking of tests that are set to be re-run

Results Review Improvements

- The option to display reportable results only has been added to the Review Results Confirmation screen, and you can now also view previous results for the selected accession(s) from the Review Results Confirmation screen
- Audit details for previous results can now be accessed from the Previous Result column of Review Results

Micro Enhancements

- The option to suppress isolates has been added to the Review Results screen



- Microbiology tests with the type General will appear in their own section of the Review Results Confirmation screen with each isolate displaying in its own column and all linked results appearing below the isolate

Workflow Management Improvements

- The Plate Reagent Lot button has been renamed **Plate Inventory**, from which Consumables and Equipment can now be added, along with Reagents, for inventory tracking

For cannabis testing laboratories:

In addition to the Workflow Management Improvements listed above:

Improved LOD/LOQ Results Translations

- This improvement will ensure better integration with third parties such as METRC and Confident Cannabis

In addition to these enhancements, there are a number of cumulative features that have been delivered throughout multiple versions of CGM LABDAQ:

- An outline of these is set out in **Appendix A**.
- **Appendix B** provides answers to Frequently Asked Questions on our recent switch from Oracle to Microsoft that was completed in March of 2020, and delivered in version 20.1

We look forward to offering this version, and continuing our commitment to supporting your operations.

Sincerely,

Carl Smith

GM, CompuGroup Medical North America Lab Division



CompuGroup Medical US
10065 Red Run Blvd., Ste. 150
Owings Mills, MD 21117

T +1 (443) 334-1830
M +1 (410) 303-3803
F +1 (443) 817-0816

carl.smith@cgm.com
<https://www.cgm.com/us>