



Dear Valued Customer,

I am pleased to announce the General Availability of CGM LABDAQ® ver. 22.1 from CompuGroup Medical.

The key improvements with this release that have been introduced since ver. 21.3 was released in 2021 include:

Gender Support Improvements

- As various states require labs and EMRs to report on additional gender data (for example where a patient's gender identity is different from their biological gender and both must be captured correctly), CGM LABDAQ has been extensively enhanced to support multiple genders.
 - CGM LABDAQ users can now define their own genders, to accurately reflect a patient gender based on their state or other specific requirements
 - This support matches Medicare requirements
 - This support is extensive throughout CGM LABDAQ, including in:
 - General set up
 - Test, facility, location, rules, and patient set up
 - Reports
 - Importing and exporting

Result Flagging Improvements for Priority Results

- Flag Icons are now provided on panels that contain out of range results to enable staff to see immediately without having to expand each panel
- Panels in Requisition Maintenance are now flagged the same as the test with the highest severity level
- Manager Review will also display flags on the panel level

Scanned Documents can now be created as PDFs

- Scanned documents can now be created as PDF's through any Attachments section within CGM LABDAQ

Control Alert Settings Enhancements

- Users in a specific user group can now see control alerts for all locations (not just their logged in location)

Performance Enhancements

- "Maximum Thread Count" has been added to the exporter to help improve the speed of each export

Administration Enhancements



- Four new housekeeping jobs have been added, including: Delete Export History rows over 180 days old, Delete Import Log rows over 180 days old, Delete Dispatch Queue Labels over 7 days old, and Delete Dispatch Queue Reports over 180 days old
- Support menu enhancements
 - The support menu has been moved from the CGM LABDAQ logo to the System Functions section
 - A link to CGM LABDAQ video tutorials has been added
 - Home page tiles can now be created for each item

General Enhancements

- Sets have been added as criteria for many rules
- Custom badges can be added to panels and tests for easy identification

Specifically for cannabis laboratories:

Web Services: Signatories

- Users can now choose one of two roles for the person signing the Certificate of Analysis
 - “Run-by User”
 - “Approved-by User”
- Users can also choose whether to use the signatory’s initials or full name

Confident Cannabis Signature URL

- Users can store the Confident Cannabis Signature URL so the user’s signature that is kept in Confident Cannabis appears on the Certificate of Analysis

Weight on Hand

- Confident Cannabis can return “Weight on Hand”
 - CGM LABDAQ can be set up to calculate a sample’s remaining weight after the amount used for a test is entered

In-Progress Orders

- “In Progress” Status will import from Confident Cannabis

In addition to these enhancements, there are a number of cumulative features that have been delivered throughout multiple versions of CGM LABDAQ. See [Appendix A](#).

[Appendix B](#) provides answers to Frequently Asked Questions regarding the recent switch from Oracle to Microsoft SQL Server that we introduced in March of 2020, with CGM LABDAQ ver. 20.1.

We look forward to offering this version and continuing our commitment to supporting your operations.



Sincerely,

Carl Smith

General Manager, Lab Division



CompuGroup Medical US
10065 Red Run Blvd., Ste. 150
Owings Mills, MD 21117

T +1 (443) 334-1830
M +1 (410) 303-3803
F +1 (443) 817-0816

carl.smith@cgm.com
<https://www.cgm.com/us>