



Dear Valued Customer,

I am pleased to announce the General Availability of CGM LABDAQ® ver. 23.1 from CompuGroup Medical.

The key improvements with this release include the following enhancements to these available software options: Lab Insights (formerly Manager Review), PCR Linking, LAB IQ, and Inventory Management.

Lab Insights with Turnaround Time Monitoring for Pending Samples

- Turnaround Time information can now be added to the Requisitions tab of Lab Insights
- This enables lab managers and directors to review on one screen whether their samples that are in pending status are predicted to meet their target Turnaround Time (as set by the lab)
- This also provides Turnaround Time information for panels that have been run that require approvals

Perform Microbiology Tests on a PCR Instrument

- CGM LABDAQ now includes the ability to easily use PCR instruments to perform microbiology testing, and identify Isolates and their Sensitivity/Resistance Based on PCR results
- CGM LABDAQ can now report isolates and medication sensitivity, after ruling out known genetic resistance markers, from PCR testing.
- This is enabled by a new PCR Linking option

LAB IQ: Profit View

- Our LAB IQ analytic module now includes a Monthly Financial Trends capability that shows the profitability of each panel by displaying reimbursement minus cost.

Improvements to Inventory Management

- Barcodes can now be scanned to add new lots in Inventory Maintenance.
- Inventory looks up by manufacturer is now available
- Email Alerts can now be set up to bring attention to low inventory

In addition to these enhancements, there are many cumulative features that have been delivered throughout multiple versions of CGM LABDAQ. See [Appendix A](#).

[Appendix B](#) provides answers to Frequently Asked Questions regarding the update from Oracle to MS-SQL Server that we introduced in March 2020, with CGM LABDAQ ver. 20.1.

We look forward to offering this version and continuing our commitment to supporting your operations.

Sincerely,

Carl Smith

General Manager, Lab Division



CompuGroup Medical US
10065 Red Run Blvd., Ste. 150
Owings Mills, MD 21117

T +1 (443) 334-1830
M +1 (410) 303-3803
F +1 (443) 817-0816

carl.smith@cgm.com
<https://www.cgm.com/us>