



Valued Clients,

I am pleased to announce the General Availability of CGM LABDAQ v24.4 from CompuGroup Medical.

The key improvements introduced to CGM LABDAQ Laboratory Information System in v24.4 are summarized below.

Auditing improvements

- Includes the ability to search and view the audit log for reporting type (to see only physician reports, for example)

Improvements to the PCR Linking option

- This software option supports the use of PCR tests for microbiology tests. New improvements to the module include:
 - Ability to suppress PCR Analyzer Results on a panel (since not all isolates tested for may be on a particular panel) and confine the results to panel-specific organisms
 - If an isolate result is rejected, all its sub-tests will automatically be rejected

Rules expansion for more robust microbiology and billing use

- Microbiology Rules improvements (also relevant to the PCR Linking option) where users can now leverage the following:
 - The new criteria of Result Flag Result Qualifier
 - The new action of Unsuppress an isolate
- Billing Rules improvements
 - Users can now leverage Result Flag and Result Qualifier

Interface Management improvements

- Users can see Work Order information when looking up Work Orders by accession in the reconciler, plus more information about the panel
- A new Error Counts Alert enables users to set an alert for import and export errors that exceed x errors per y minutes for prompt troubleshooting

In addition to these enhancements, there are many cumulative features that have been delivered throughout multiple versions of CGM LABDAQ. See [Appendix A](#).

[Appendix B](#) provides answers to Frequently Asked Questions regarding the update from Oracle to MS-SQL Server that we introduced in March 2020, with CGM LABDAQ v20.1.

We look forward to offering this version and continuing our commitment to supporting your operations.

Sincerely,



Carl Smith

COO, CompuGroup Medical



CompuGroup Medical US
10065 Red Run Blvd., Ste. 150
Owings Mills, MD 21117

T +1 (443) 334-1830
M +1 (410) 303-3803
F +1 (443) 817-0816

carl.smith@cgm.com
<https://www.cgm.com/us>