

Product: DAQbilling**Release Date:** January 2015**Customer Support:** Contact the Help Desk at (800) 359-0911 option 5 with any questions.

Introduction

This document contains information on gaining access to CompuGroup Medical DAQbilling software. Please forward to your IT department. To access the software, you need a computer with internet access and the following specifications.

Minimum Specifications

- 900 MHz CPU
- 256 MB RAM
- 500 MB Free Disk Space
- Windows 7, Windows 8/8.1 Pro, Vista or above
- Stable high speed internet connection (DSL or broadband cable)
- Monitor – requires a screen resolution of at least 800x600 with 16-bit color

Recommended Specifications

- 2.0 GHz CPU
- 512 MB RAM
- 1 GB Free Disk Space
- Stable high speed internet connection (DSL or broadband cable)
- Windows 7, Windows 8/8.1 Pro, Vista or above
- Monitor – requires a screen resolution of at least 800x600 with 32-bit color
- Color Laser Printer
- Scanner – must be twain compliant and 32-bit compatible
- **For sites with firewalls, the PC needs to access to these IP addresses: 74.217.186.11, 74.217.186.12, 74.217.186.13, 74.217.186.14**
- **TCP Ports 9020 – 9040 need to be open to all IP addresses**
- **Interface connections need access to the IP address 74.217.186.10, and Ports 9001-9002 need to be open**